**SE Manual BigJack IR-6SP**

**GB Manual BigJack IR-6SP**

**DE Manual BigJack IR-6SP**

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# GB Manual BigJack IR-6SP

## *Introduction*

BigJack is a large and robust remote control unit to which external control switches can be connected. It can be programmed, which means that it is possible to "teach" it IR codes from other remote control units, such as those used for TVs, videos, stereo units and toys. It is suitable for use by people who have cognitive and motor difficulties, and who require large and clear control buttons.

## *General Introduction to BigJack*

There are certain things you should know before starting to use BigJack.

- In order to enter the programming mode, it is necessary to press the programming button (located underneath the unit) with, for example, a pen while **at the same time** pressing a function button. Different function buttons are used depending on what is to be programmed.
- You can always press the programming button if you press the wrong button or if something goes wrong with the programming. The programming button functions as a reset button.
- It is important not to wait too long between the different steps of the programming procedures. BigJack automatically exits from programming mode after a certain period in order to save power.

The figure shows different names and functions that are described in the text.

Fig 1



Fig 2

Certain basic settings have the following factory defaults when BigJack is delivered:

- The sound volume is set to 3.
- The response delay is set to 1, meaning no delay.
- GewaLink channels 16-21 are programmed to the buttons 1-6.

## Description of Symbols

The following symbols will be used in the text:

- **P** = programming button
- **1** .... **6** = function buttons 1-6

## Inserting the Batteries

Unscrew the screws as shown in Figure 3 using a Philips screwdriver. Remove the upper part and insert four 1.5 V alkaline batteries of type LR6. Insert the batteries with the correct polarity as shown on the bottom of the battery compartment. It is important to use **alkaline** batteries that are protected from leakage and that have a long lifetime. Finish inserting the new batteries by **pressing** the programming button.

Fig 3



Fig 4

## Programming

### Teaching IR Codes from Another Remote Control Unit


In order for BigJack to control a TV set, video, stereo unit, DVD, toys, lamps, etc., it is necessary that these appliances have a remote control that emits IR. Most remote control units on the market use infrared light (IR), as does BigJack. A TV set, for example, always has a remote control as standard, which is used to control various functions. It is this remote control unit that is to be used when teaching BigJack the IR channels. Other appliances such as lamps, hall telephones, telephones and automatic door openers do not have such a remote control unit. Gewa AB supplies solutions that make it possible to control such appliances. Contact Gewa AB or a representative of Gewa in such cases.

Proceed as follows to teach BigJack IR codes from another remote control unit:

1. Check that the batteries in BigJack and in the remote control unit are fully charged. Shield the units from intense light such as direct sunlight during the teaching procedure.
2. Place BigJack on a table with the remote control unit pointing at the IR sensor on BigJack, as shown in Figure 5.
3. Press **P** and **1** at the same time.

In order to check that BigJack and the remote control unit are positioned correctly, press and hold down any button on the remote control unit. Move the remote control unit left and right until the green indicator lamp goes out. Then move the remote control unit backwards until the indicator lamp goes out (1 metre for certain units). Place the remote control unit between these limits. Certain types of remote control unit give only a short flash for each press. If this is the case, press the button on the remote control unit repeatedly.

Fig 5

4. Press and release the button on BigJack on which the IR channel is to be placed. The red and green indicator lamp now lights with a steady light.
5. Press and hold down the button on the remote control unit that is to be taught to BigJack. Release it when the red indicator lamp goes out. This may take as long as 10 seconds. If the red indicator lamp returns to a red light, press again one or several times until both red and green indicator lamps go out. If the indicator lamp flashes red and a sound signal is heard ♪♪♪♪ , the teaching process has failed. Repeat the process from step 4.
6. Conclude the process by pressing button **6**.
7. If the teaching process has succeeded for the first button, continue with the next button as described above in points 3-6 until all IR channels are programmed.

**Hint!** BigJack can learn IR channels from most remote control units on the market, but there are a few exceptions. It may affect certain functions or it may affect all of the functions on a remote control unit. It is possible in these cases to use another teaching mode. Replace step 3 above by the following procedure: press in **P** and **1** at the same time and hold them depressed for about 2 seconds until the green indicator lamp flashes and then goes out. The unit is now in a different learning mode. Continue with the procedure described above.

### **Quick guide**

Teaching: **P** + **1**, choose button on BigJack, press the button of the remote control unit, conclude with **6**.

## Reset to Factory Defaults

It is possible to reset all of BigJack's buttons and settings to the factory defaults, shown below:

- Sound volume reset to value 3.
- Delay time reset to value 1, no delay.
- GewaLink channels 16-21 programmed to the function buttons 1-6.

Proceed as follows to reset to factory defaults:

1. Press **P** and **3**.
2. Press the following buttons in this order **5**, **1**, **2**, **3**, **6**. The indicator lamp on button 1 flashes green to show that the process is complete.

### **Quick guide**

Reset to factory defaults:

**P** + **3**, **5**, **1**, **2**, **3**, **6**

## Setting the Sound Volume

BigJack has a sound signal that is activated when a button is pressed. The volume of this signal can be regulated in six steps (**1**, **2**, **3**, **4**, **5**, **6**).

**1**= no sound signal and **6** = loudest sound signal.

Proceed as follows:

1. Press **P** and **4** at the same time.
2. Press a button from **1** - **6** to select the sound volume. The indicator lamp on button 1 flashes green to show that the process is complete.

### **Quick guide**

Set the sound volume: **P** + **4**, **1** - **6**.

**1**= no sound signal and **6** = loudest sound signal.

## Setting delay time

It is possible to set the delay time for BigJack, the time that a button must be held down before a response is obtained. The delay time can be set to one of six values (1, 2, 3, 4, 5, 6).

1 = no delay time and 6 = maximal delay time (0.5 sec).

Proceed as follows:

1. Press **P** and 2 at the same time.
2. Press a button from 1 - 6 in order to select the time of the delay time. The indicator lamp on button 1 flashes green to show that the process is complete.

### Quick guide

Set the delay time: **P** + 2, 1 - 6.

1 = no delay time and 6 = maximal delay time.

## Learning sequences (Macro)

BigJack is able to send sequences of signals. Thus it is possible for BigJack to send a series of button presses following the pressing of **one button**. One example of this is switching off the sound of the TV when answering the telephone. BigJack can remember the length of time that a button is to transmit.

It is always possible to interrupt the transmission of a sequence by pressing any other button. The maximum sequence that can be learnt lasts for 25.2 seconds and contains up to six presses on each button.

NOTE! Buttons 5 and 6 cannot store sequences, and the button on which a sequence is stored cannot be included in the sequence.

1. Press **P** and 5 at the same time.
2. Press 1.
3. Press the buttons that are to be included in the sequence. Hold them pressed for the length of time that they are to transmit when the sequence is replayed.
4. Press **P** and 5 at the same time.
5. Press 2.

6. Select the button on which the sequence is to be stored (do not select 5 or 6). The indicator lamp on button 1 flashes green to show that the process is complete.

### **Quick guide**

Learning sequences: **P** + **5**, **1**, press the sequence, **P** + **5**, **2**, select the button for the sequence.

## **External Control Switches**

It is possible to connect external control switches to BigJack, which influence buttons 1-6. The control contacts are connected to the inputs on the front using 3.5 mm telephone plugs.

## **Cleaning**

BigJack may be cleaned with a rag that has been immersed in water with a small amount of washing up liquid and then **wrung dry extremely thoroughly**. BigJack may be dried with a dry rag after cleaning.

## **Summary of Quick guide**

### **Learning channels**

Learning: **P** + **1**, select button on BigJack, press on the remote control unit, conclude with **6**.

### **Setting delay time**

**P** + **2**, **1** - **6**.

**1** = no delay time and **6** = maximal delay time.

Factory default = no delay time.

### **Reset to factory defaults**

**P** + **3**, **5**, **1**, **2**, **3**, **6**

### **Set the sound volume level**

**P** + **4**, **1** - **6**.

**1** = no sound signal and **6** = loudest sound signal.

Factory default = 3.

### **Learning sequences**

**P** + **5**, **1**, press the sequence, **P** + **5**, **2**, select button for the sequence.

## ***Accessories***

### **Accessories included with BigJack**

<b>Art. no.</b>	<b>Name</b>	<b>Description</b>	<b>Quantity</b>
-----	Note sheet BigJack	A4 sheet for recording programmed channels	1

### **Other accessories for BigJack**

<b>Art. no.</b>	<b>Name</b>	<b>Description</b>
8362	Battery LR6	1.5 V Alkaline

## ***Specifications***

Current drawn during stand-by:	approx. 10 $\mu$ A
Current drawn transmitting:	approx. 35 mA (GewaLink channel)
Battery life:	approx. 2 years (200 seconds transmission per day on GewaLink channels)
Programming locations:	6
Operating force for buttons:	approx. 1.5 N (150 gram)
Range:	20-30 m (GewaLink channels)
Dimensions (L x W x H):	140 x 100 x 38 mm
Weight with batteries:	280 gram
Temperature range (not valid for batteries):	-20°- +45°C
Contacts for connecting control contacts:	3.5 mm

## ***Compliance with the Medical Device Directive***

IR-6SP BigJack is marketed as a technical aid for persons with physical disease and it satisfies the requirements of the Medical Device Directive 93/42/EEC. These include requirements for EMC and for electrical safety.

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